

Photopatterning of Methacrylated Ligands for Spatial Control of SynNotch Receptors

Alexander Do, Stephaine Do, Megan McCain

Dept of Biomedical Engineering, Bridge Institute, University of Southern California, Los Angeles, CA, USA



Bridge UnderGrad Science (BUGS) Summer Research Program

Abstract

- Engineered human tissue doesn't allow spatial control of cells in multi-cell tissues and results in poor representation of true physiology.
- I spatially controlled synNotch activated cells by UV-activated crosslinking of chemically modified GFP onto Gelma.
- The crosslinking resulted in a high density of GFP in the region of UV exposure and cellular activation occurred more often in the area.

Figure 1





Figure 2



Methods

Figure 3A




Figure 3B





Figure 3C



Figure 3D



- Gelma is added to the exposed section of the coverslip.
- PDMS is added on top to stamp down gelma
- GFP and photoinitiator is mixed in and placed on the gel.

- Photomasks are made from laser engraved taped glass.
- Coverslips are placed into 3D mold with photomask layed on top.

- Device from Figure 3B is placed into light box.
- Time and intensity determined by Figure 5.

- GFP is crosslinked to gelma in UV exposed region(see Figure 2).
- Fibroblasts are seeded on top with activation occurring when in contact with the GFP.

Results

Figure 4A




Figure 4B




Figure 4C




Figure 5




Fig 4. Activated SynNotch cells on photopatterned GFP Fig. 4A. GFP Fig. 4B. mCherry Fig. 4C. Brightfield

Conclusion

- This process will eventually lead to controlled cell differentiation to represent tissues of multiple cell types like neuromuscular junctions.
- The modification can be used on any protein and future studies can go beyond synNotch and cellular differentiation.
- The impact of more physiological accurate engineered tissues is faster drug production, less dependency on animal models, and eventually lab grown organs.

CONTACT US

bridge.usc.edu/bugs,
doalexander700@gmail.com